



- ① = Option: filtering system (exception: above the fume hood), ideally a separate filter cabinet containing 3 to 4 filtering stages is placed next to the fume hood
- ② = Option: interior lining in nickel chromium steel or polypropylene
- ③ = Front sash made from Macrolon/acrylic glass in one part, Option: Openings for reaching inside
- ④ = Service fitting points
- ⑤ = Aluminium channel with electrical fixtures and remote control fittings
- ⑥ = Option: heavy-load support frame
- ⑦ = Option: screen against gamma and beta rays