


Fume cupboard systems | 1.5 DELTAguard filter fume cupboard



Variation: Overall depth 1070 mm, Useable depth working surface 880 mm instead of 720 mm

- | | |
|---|---|
| 01 Filter | 07 Aluminium duct with electrical service outlets and front control valves |
| 02 Option: sash with 3/4 cross slides | 08 Option: Under bench cabinet or safety store cabinet |
| 03 Fume cupboard air functional monitoring display | 09 Option: Drip cup |
| 04 2 rows tripod rod holder | 10 Lighting |
| 05 Option: Internal sockets, switched outside | |
| 06 Service fitting points | |

- Available as one or two phase filter
- 2 phase filter from frame size 1500 (asymmetrical sequence)