


## Table systems | 3.4 Auxiliary frame



Auxiliary frame

With the auxiliary frame, a metal frame is the supporting element of the worktop. It is covered and has a 100 mm base. The auxiliary frame can be used as a stand-alone unit (wall bench and island bench), in combination with a standing bridge and with all table systems.

- Standard widths: 600 / 900 / 1200 / 1500 / 1800 mm
- Standard heights: 750 / 900 mm
- Standard depths: 160 mm

**01** Option: drip cup

**02** Option: Media