


Fume cupboard systems | 1.9 DELTAguard Ex-fume cupboard



Variation: Total depth 1070 mm, useable depth 870 mm instead of 710 mm

- 01** Upper part glazed fixed screen
- 02** Sash with fixed glass
- 03** Option: (depending on EX-zone)
Exhaust functional monitoring display EX model
- 04** 3 series tripod rod holders from PPsEL
- 05** Option: Electrical service outlets in EX-protected installation housing
- 06** Service fitting points in explosion-proof design
- 07** Option: Aluminium duct with front control valves in explosion-proof design
- 08** Option: With base cupboard or safety cabinet in explosion-proof design
- 09** Option: Drip cup
- 10** Option: Lighting explosion-proof design