


## Table systems | 3.2 Series 4FU



## Series 4FU

In the Series 4FU a metal 4-legged frame is the supporting element of the worktop. This table can be used with a suspended under bench cabinet, and with an inserted under bench cabinet with 100 mm plinth, both as stand-alone, as well as with mobile substructures with open base (100 mm). The series can be combined with all media/table systems.

- Standard widths: 600 / 900 / 1200 / 1500 / 1800 mm
- Standard heights: 750 / 900 mm (seated / standing working heights)
- Standard depths: 600 / 750 mm